






INNOVATIVE LEADERS IN AIRWAY MANAGEMENT

THE guardianCPV LARYNGEAL MASK... THE ULTIMATE SEAL





INNOVATIVE LEADERS IN AIRWAY MANAGEMENT

THE guardianCPV LARYNGEAL MASK... THE ULTIMATE SEAL



- Single use
- Silicone
- Cuff Pilot Valve
- High seal pressure
- Flexible
- Gastric channel
- Built-in bite block
- Scavenge Port



SELECTION GUIDELINES

SIZE	3	4	5
PRODUCT CODE:	LMGCP03	LMGCP04	LMGCP05
PATIENT SELECTION	30-50kg	50-70kg	70-100kg
CUFF PILOT RECOMMENDATION	Green Zone Pressure Range = 40-60cm H ₂ O for all sizes		

Available for release Q1-2011



Distributed by
Ultimate Medical Pty Ltd
484 Swan Street, Richmond
Victoria 3121 Australia

t: +61 3 9421 0500
f: +61 3 9421 0555
www.ultimatemedical.com.au
www.cuffpilot.com