




# The Ultimate Flexy™ Laryngeal Mask with Cuff Pilot™



## FLEXY SINGLE USE

- Cuff Pilot™ pressure indicator
- Silicone Cuff for smooth insertion
- Flexible wire reinforced tube, useful in head and neck, oral and nasopharyngeal procedures
- A new mask for every patient
- Traditional sizes
- Latex free

## SELECTION GUIDELINES

SIZE	1	1.5	2	2.5	3	4	5	6
PRODUCT CODE	LMFLCP01	LMFLCP1.5	LMFLCP02	LMFLCP2.5	LMFLCP03	LMFLCP04	LMFLCP05	LMFLCP06
PATIENT SELECTION	NEONATES 0kg-5kg	INFANTS 5kg-10kg	INFANTS 10kg-20kg	CHILDREN 20kg-30kg	CHILDREN 30kg-50kg	ADULT 50kg-70kg	ADULT 70kg-100kg	ADULT 100kg+
CUFF PILOT™ RECOMMENDATION	Green Zone = P ≤ 60cm H <sub>2</sub> O (for all sizes)							

The Cuff Pilot™ is a guide to the intracuff pressure. Physicians should always use their best clinical judgement in regard to the appropriate level of inflation of the airway.

Distributed by Ultimate Medical Pty. Ltd.  
484 Swan Street, Richmond Victoria 3121 Australia  
t: +61 3 9421 0500  
f: +61 3 9421 0555

[www.ultimatemedical.com.au](http://www.ultimatemedical.com.au)  
[www.cuffpilot.com](http://www.cuffpilot.com)

